

## Experimental characterisation of a vacuum arc thruster using Impedans' Semion Pulsed DC retarding field energy analyser

### INTRODUCTION

Nanosatellites, also called CubeSats, are increasingly being utilized as teaching tools for technological demonstrations and for Earth observation and scientific applications. Because of their small size, limited available power and restricted mass budget, the current generation of CubeSats are not equipped with a means of propulsion. In cases where a conventional propulsion system is installed, it takes a significant share of the available spacecraft volume. Nevertheless, in order to further develop the use of CubeSats, propulsion is a major requirement.

Vacuum Arc Thrusters (VATs) are uniquely suited to such applications because of their simple, scalable, and robust construction. VATs do not require cathode neutralizers and use compact solid metal propellant. However, ablative pulsed thrusters, such as VATs are particularly prone to damage due to their high current operation and are typically limited by the number of pulses they can produce. A long duration type of VAT known as the inline screw-feeding vacuum-arc-thruster (ISF-VAT) addresses the destructive effect of the vacuum arc discharge on thruster components, namely, the electrodes and insulator.

### EXPERIMENT

The ISF-VAT concept, see figure 1, is built around a central cathode rod that sits within a concentric insulator tube. A second electrode, positioned at the




Figure 1. Schematic diagram of the ISF-VAT.

outer edge of the insulator, functions both as the anode of the dc circuit and as the exit plane of the thruster. The cathode is connected to a metallic headless screw, which allows the cathode to be advanced at a precise rate inside the insulator, in a helical motion, which compensates from the both the radial and azimuthal cathode erosion patterns allowing for constant thruster geometry throughout the lifetime of the thruster. The thruster plume was characterized using a commercial retarding field energy analyser (RFEA) ([Semion](#), [Impedans](#)) using a single sensor. The measurements are performed sweeping the RFEA around the thruster exit in an arc, at a distance of 150 mm in the plane of the thruster. The RFEA electronics signal sample duration is 1  $\mu$ s.

### RESULTS

The discharge current peak and duration is influenced by the quality of the cathode layer deposition on the insulator. In figure 2, three separate arc discharge currents are shown, representing the lowest, typical and the strongest individual pulse energies




Figure 2. Temporal evolution of the Normalised discharge current and total ion current.

recorded during the entire (12 hr) run. The discharge current typically lasts 200  $\mu\text{s}$ . The total ion current follows a similar decay rate as the typical discharge current. This result permits the calculation of the ion fraction of the discharge current as  $f_i = I_{\text{ion}}(t)/I_d(t)$ . At  $t_p = 20 \mu\text{s}$ , the absolute measured values are  $I_d(t_p - 10\mu\text{s}) = 11.5 \text{ A}$  and  $I_{\text{ion}}(t_p) = 0.57 \text{ A}$ , respectively and then  $f_i = 0.05$ , which is in the range given in the vacuum-arc literature of 0.05–0.1.

The ISF-VAT plume was characterized using the movable RFEA sensor. The data were collected over many thousands of firing pulses. The measurements were performed with the RFEA sensor located at the azimuthal angles of  $0^\circ$  to  $75^\circ$  in  $15^\circ$  increments and between  $80^\circ$  to  $95^\circ$  in  $5^\circ$  increments, where  $0^\circ$  corresponds to the thruster-axis direction. The RFEA measurements were synchronized with thruster ignition and recorded from 0 to 200  $\mu\text{s}$  in the intervals of 20  $\mu\text{s}$ . For each measurement condition (angular position relative to the thruster axis, RFEA grid voltage and delay time), an individual firing pulse was generated.

The time-resolved spatial distributions of the ion current density are shown in figure 3. As expected, the ion current decreases with time. The half-angle beam divergence, defined as 90% of the ion current, is  $\theta \approx 60^\circ$  in the first 20  $\mu\text{s}$  and increases with time. Fig 4 shows the IEDFs obtained with the RFEA with a corrected grid voltage. The IEDFs were obtained by fitting the  $I(V)$  data from the RFEA with an arctangent function and then taking the derivative of this fitted analytical function. From these results, the peak ion population, in all directions, corresponds to a beam




Figure 3. Ion-current density as a function of the azimuthal angle where  $0^\circ$  represents the direction the thruster is pointed.




Figure 4. IEDF for different azimuthal positions ranging from  $-5^\circ$  to  $90^\circ$ .

kinetic energy of ( $E_{\text{ion}} \approx 65 \text{ eV}$ ). Using the ion energy the most probable ion velocity can be calculated as follows:

$$v_i = \sqrt{(2(E_{\text{ion}})/M_i)} \approx 16 \text{ km/s,}$$

Where the mass of the ion is  $M_i = 7.95 \times 10^{-26} \text{ kg}$ . The time of flight from the thruster to RFEA probe is  $\approx 10 \mu\text{s}$ .

## CONCLUSION

An ISF-VAT was operated continuously for 12 hrs permitting the characterisation of the thruster over its lifetime. The ISF-VAT parameters were characterised using the RFEA measurements of the plume.

## REFERENCES

\*Kronhaus, I. et al, "Experimental Characterisation of the Inline-Screw-Feeding Vacuum-Arc-Thruster Operation"

IEEE Transaction on Plasma Science.

Doi: 10.1109/TPS.2017.2776839